

Bed Alarm & 6 Month Sensor Pad




The Avalon Bed Alarm System is a device that is used to monitor a patient in a wheelchair or bed to prevent unintentional falls. The Bed Alarm System consists of two main components, an Alarm and Sensor pad. The Sensor pad is pressure sensitive, so as the patient rises off the pad, the alarm is set off. The Sensor pad is made-up off tiny magnetic sensors throughout the pad that are attached to a seven foot (phone style) cord. The tip of the cord has a phone jack connector that plugs in to the bottom of the alarm.

The alarm monitor works in conjunction with the sensor pad. The Alarm has an audible, beeping sound that alerts caregivers that the patient is leaving there bed or wheelchair. The alarm has a High/Low setting as well as an on/off switch. If you are in a Nursing Home or Hospital with a Nurse call system, the Avalon Bed Alarm and Sensor pad system also comes with a nurse call cord and two-in-one adaptor. The portable patient alarm can function on AC/DC Power (batteries or Plug).



This alarm System is also known as: Bed Alarm System, Patient bed alarms, bed alarms for seniors, elderly bed alarm, bed alarm system, bed alarm for elderly, and fall prevention for seniors.



Fall Prevention is a major problem in the senior popular today. There are many factors that can effect person coordination such as medication, illness, eyesight, hearing, coordination, balance, and physical strength. As a person grows older, your bones become more brittle, and your body heals and fights off disease at a slower rate. Falls are the leading cause of injury and accidental death for older adults, so with a few preventive measures like a bed alarm system, inconjunction with close monitoring these falls can be prevented. Sometimes it's hard to monitor a person 24 hours a day, and that is where the Avalon Patient Fall Prevention Alarm System comes in to assist you!

For more State-of-the-art Fall Prevention products for seniors, please visit our website for the Lowest pricing on the web!

Mountainside Medical Equipment, Inc.

9262 Old River Road

Marcy New York 13403

Phone: 315-768-3030

Fax: 315-768-6111

www.Mountainside-Medical.com