

SECTION 2—MODEL NUMBERS

MODEL NOS. STD.	DESCRIPTIONS
9042	Polyester Sling
9043	Polyester Sling With Commode Opening
9046	Polyester Mesh Sling
9047	Polyester Mesh Sling With Commode Opening
9070	Strap Kit
9071	Chain Kit

Note: Slings listed above require a strap or chain kit for hook-up.

CHAIN-FREE ALL FABRIC SLINGS	
*R110 - R117	Full Body Slings With or Without Commode Openings
*R101 - R102	Divided Leg Slings
*R121	Toileting Sling
**R140 AND R141	Heavy Duty With or Without Commode Openings
***R130, R131, R134, R136	Standing and Transport Slings
<p>*Note: Can ONLY be used with Invacare Models 9700, 9701, 9805, 9806, 9807 and Reliant Lifts on the large hooks of the swivel bar or on converted model nos. 9801, 9802 and 9803 (9800A Boom Retrofit Kit). The small hooks are to be used for chain/strap kits only.</p> <p>**Note: Can ONLY be used with Invacare Model Reliant 450 and 600 Lifts.</p> <p>***Note: Can ONLY be used with Invacare Model Reliant 350 and 440 Lifts.</p>	

WARNING

DO NOT exceed maximum weight limitation of the patient lift. The weight capacity is limited to the lowest rated capacity of any one of the components in use (e.g. patient lift, sling, or scale). The patient's weight **MUST** not exceed the lowest rated capacity of any component.