

The FILAC™ 3000 Electronic Thermometers feature completely separate color-coded isolation chambers with individual probes, probe wells and probe cover holders. This patented system provides the benefits of electronic thermometry while at the same time reducing the risk of cross-contamination.

PRODUCT SPECIFICATIONS

Accuracy

- Quick Mode $\pm 0.3C$ ($\pm 0.6F$)
- Standard Mode $\pm 0.1C$ ($\pm 0.2F$)
- Direct Mode $\pm 0.1C$ ($\pm 0.2F$)

Measurement Response Time

- Quick Mode: <4 sec Oral, <15 sec Axillary/Rectal
- Standard Mode: <10 sec Oral; <15 sec Axillary/Rectal
- Direct Mode: <60 sec Oral, Axillary, Rectal
- Cold Mode: <15 sec Oral; <20 sec Axillary/Rectal

Ambient Operating Range

- 10C - 40C (50F - 104F) at 10% - 95% RH

Temperature Range

- 30.0C - 43.0C (86.0F - 109F)

Storage Temperature Range

- -25C - 55C (-13F - 131F) at 95% RH (non-condensing)

Battery Life

- 6000 measurements on 4 'AA' Alkaline 1.5V batteries

Anti-Theft System (AD Only)

- Electronic system (return to base message)
- Mountable base with key lock
- Mobile rolling stand base with key lock

Standards

- UL60601-1 Compliant
- CAN/CSA C22.2No601.1 Compliant
- EN 12470-3 (In standard predict mode)
- ASTM 1112 (In standard predict mode)

ORAL/AXILLARY ISOLATION CHAMBER



RECTAL ISOLATION CHAMBER

Complete Units

Code	Description	Ship Case
504000	FILAC 3000 EZ Oral/Axillary Complete	1 Each
504005	FILAC 3000 EZ Rectal Complete	1 Each
505000	FILAC 3000 AD Oral/Axillary Complete	1 Each
505005	FILAC 3000 AD Rectal Complete	1 Each

Probes and Chambers

Code	Description	Ship Case
500026	FILAC 3000 Oral Probe-4ft	1 Each
500027	FILAC 3000 Oral Probe-9ft	1 Each
500028	FILAC 3000 Oral Isolation Chamber	1 Each
500036	FILAC 3000 Rectal Probe-4ft	1 Each
500037	FILAC 3000 Rectal Probe-9ft	1 Each
500038	FILAC 3000 Rectal Isolation Chamber	1 Each

Probe Covers

Code	Description	Ship Case
500500	FILAC 3000 Probe Covers	500 Each/Case
502000	FILAC 3000 Probe Covers	2,000 Each/Case

Accessories

Code	Description	Ship Case
500046	FILAC 3000 Wall Mounting Base with Lock	1 Each
500047	FILAC 3000 Mounting Base for Rolling Stand	1 Each
500048	FILAC 3000 Mounting Base with Pole Clamp	1 Each
500099	FILAC 3000 Calibration Plug	1 Each
303059	Rolling Stand	1 Each
8884813711	Rolling Stand Basket	1 Each

COVIDIEN, COVIDIEN with Logo and ™ marked brands are trademarks of Covidien AG or its affiliate.
©2007 Covidien AG or its affiliate. All rights reserved.



H-6122-10M-01/08

15 HAMPSHIRE STREET
MANSFIELD, MA
02048

1-800-962-9888
508-261-8000

WWW.COVIDIEN.COM