



COVIDIEN

positive results for life™

Filac™ 3000 EZ Electronic Thermometer



Accurate:

- +/- 0.3° C Quick Mode
- +/- 0.1° C Standard Mode
- +/- 0.1° C Direct Mode

Fastest Predictive Readings:

- < 4 seconds Oral
- 10-15 seconds Axillary/Rectal
- < 60 seconds Direct Mode

Easy to Use:

- Intuitive icon based user interface
- Easy to read LCD display with backlight
- Durable, robust design

Supports Infection Control:

- Interchangeable, color coded isolation chambers
- Self contained isolation chamber design
- Single use durable probe covers minimize cross contamination



Accuracy, At Your Pace. *Fast*

PRODUCT SPECIFICATIONS

ACCURACY

- Quick Mode: (Oral) +/- 0.3° C (+/-0.6° F)
- Standard Mode: (Axillary/Rectal) +/- 0.1° C (+/-0.2° F)
- Direct Mode: +/- 0.1° C (+/-0.2° F)

MEASUREMENT RESPONSE TIME

- Quick Mode: <4 second Oral
- Standard Mode: 10-15 second Axillary Rectal
- Direct Mode: <60 second Oral, Axillary, Rectal

AMBIENT OPERATING RANGE

- 10° C - 40° C (50°F - 104°F) at 10%-95% RH

TEMPERATURE RANGE

- 30° C - 43.0° C (86.0°F - 109°F)

STORAGE TEMPERATURE RANGE

- -25° C - 55° C (-13°F - 131°F) at 95% RH (non-condensing)

BATTERY LIFE

- 6,000 Temperatures (4AA Alkaline 1.5V batteries)

Order Information

Thermometer

Code	Description	Ship Case
504000	Filac 3000 EZ Oral/Axillary Complete	1
504005	Filac 3000 EZ Rectal Complete	1
504027	Filac 3000 EZ Oral/Axillary Complete 9ft	1
504037	Filac 3000 EZ Rectal Complete 9ft	1
505000	Filac 3000 ADOral/Axillary Complete	1
505005	Filac 3000 AD Rectal Complete	1
505027	Filac 3000 EZ Oral/Axillary Complete 9ft	1
505037	Filac 3000 AD Rectal Complete 9ft	1

Probes and Chambers

500026	Filac 3000 Oral Probe 4ft	1
500027	Filac 3000 Oral Probe 9ft	1
500028	Filac 3000 Oral Insolation Chamber	1
500036	Filac 3000 Rectal Probe 4ft	1
500037	Filac 3000 Rectal Probe 9ft	1
500038	Filac 3000 Rectal Insolation Chamber	1
500500	Filac 3000 Probe Covers	500 Each/Case
502000	Filac 3000 Probe Covers	2000 Each/Case

Accessories

500046	Filac 3000 Wall Mounting Base with Lock	1
500047	Filac 3000 Wall Mounting Base for Rolling Stand 1	1
500048	Filac 3000 Mounting Base with Pole Clamp	1
500099	Filac 3000 Calibration Plug	1
303059	Filac 3000 Rolling Stand	1
500099	Filac 3000 Rolling Stand Basket	1

For additional information contact Customer Service at 1-800-962-9888



COVIDIEN, COVIDIEN with Logo and TM marked brands are trademarks of Covidien AG or an affiliate. ©2008 Covidien AG or an affiliate. All rights reserved.

H6282 10M 07/08