



**COVIDIEN**

*positive results for life™*

## Genius™ 2 Infrared Tympanic Electronic Thermometer

### FEATURES

- Durable tip based sensor technology
- Tapered probe design
- Rapid response time
- Real time temperatures
- Programmable display
- Single use disposable probe covers

- The Genius 2 infrared tympanic electronic thermometer is the gold standard in thermometry
- Utilizing our patented Peak Select System, the Genius 2 thermometer records, analyzes and selects the highest and most accurate temperature in less than two seconds
- The Genius 2 thermometer optimizes patient care by delivering fast, accurate temperatures and supporting infection control initiatives





# Smarter Than Ever

## PRODUCT SPECIFICATIONS

### ACCURACY

- +/- 0.1° C (36.7 - 38.9)

### MEASUREMENT RESPONSE TIME

- 1 - 2 Seconds

### AMBIENT OPERATING RANGE

- 16° C - 33° C (60.8F - 91.4F)

### TEMPERATURE RANGE

- 33° C - 42.0° C (91.4F - 107.6F)

### STORAGE TEMPERATURE RANGE

- -25° C - 55° C

### BATTERY LIFE

- 15,000 Temperatures (3AAA batteries)

### ANTI-THEFT SYSTEM

- Wall mount bracket/Genius™ 2 mobile cart

### STANDARDS

- ASTM and CEN compliant

### WARRANTY

- 1 year full manufacturer warranty

## Genius 2 Thermometer

- Peak Select technology records and analyzes up to 100 temperature measurements and selects the most accurate reading in 1 to 2 seconds.
- Re-designed tapered probe tip facilitates proper ear canal alignment minimizing technique related errors and increasing patient comfort.
- Rugged tip based sensor construction improves durability and accuracy, minimizing the need for recalibration.
- No touch probe cover application and ejection reduce potential cross contamination.
- Programmable display shows oral, rectal, axillary and core or an absolute temperature reading in ear mode.

## Order Information

Code	Description	Ship Case
303000	Genius 2 Thermometer with Base	1
303030	Genius 2 Probe Cover	2112 each/case
303058	Genius 2 Locking Mount Wall	1
303059	Genius 2 Cart with Locking Mount	1
303060	Genius 2 Cart Conversion Kit	1
303099	Genius 2 Field Calibration Checker	1



COVIDIEN, COVIDIEN with Logo and ™ marked brands are trademarks of Covidien AG or an affiliate. ©2008 Covidien AG or an affiliate. All rights reserved.

H6283 10M 07/08