















Respond to the wound healing environment with ALLEVYN[®]


	Light Exudate	Moderate to High Exudate	High to Extra-High Exudate
Non-Adhesive	 <p>ALLEVYN Non-Adhesive For general use on wounds requiring thick, padded protection</p>	 <p>ALLEVYN Non-Adhesive For general use on wounds requiring thick, padded protection</p>  <p>ALLEVYN Heel For heels, crevices, and other difficult-to-fit or protruding areas requiring form-fitting, padded protection</p>  <p>ALLEVYN Tracheostomy For use around tracheostomy and other tubes requiring form-fitting healing</p>	 <p>ALLEVYN Cavity For deep cavity wounds that require packing; maintains shape for easier removal</p>
Gentle Adhesive	 <p>ALLEVYN Thin For fragile skin, tears, abrasions, and other wounds where a highly flexible dressing is required</p>	 <p>ALLEVYN Compression For a variety of wound types where fast exudate control is required</p>	 <p>ALLEVYN Plus Cavity For cavity, donor sites, or other wounds where persistent absorption and border conforming coverage are required</p>
Secure Adhesive	 <p>ALLEVYN Island For wounds on flat surfaces requiring a waterproof dressing and flexible protection of the peri-wound</p>	 <p>ALLEVYN Adhesive For general use on wounds where a waterproof dressing with adhesive borders to stay in place is required</p>  <p>ALLEVYN Sacrum For wounds in the sacral area requiring a waterproof covering and bacteria-proof protection with adhesive borders to stay in place</p>	 <p>ALLEVYN Plus Sacrum For wounds in the sacral area with high to extra-high levels of exudate where a waterproof covering and bacteria protection is required with adhesive borders to stay in place</p>  <p>ALLEVYN Plus Adhesive For wounds with high to extra-high levels of exudate where a waterproof covering is required with adhesive borders to stay in place</p>


● For the most common types of wounds



Ordering information

Non-Adhesive


ALLEVYN Non-Adhesive					
	Code	Description	Pcs/Pkg	Pkgs/Unit	Order Unit
	66027643	2" x 2"	10	6	CA
	66027637	4" x 4"	10	7	CA
	66020093	6" x 6"	10	3	CA
	66027638	8" x 8"	10	2	CA


ALLEVYN Heel					
	Code	Description	Pcs/Pkg	Pkgs/Unit	Order Unit
	66007630	Heel	5	6	CA


ALLEVYN Tracheostomy					
	Code	Description	Pcs/Pkg	Pkgs/Unit	Order Unit
	66027640	3 1/2" x 3 1/2"	10	8	CA

ALLEVYN Cavity					
	Code	Description	Pcs/Pkg	Pkgs/Unit	Order Unit
	66027326	2" Circular (3g)	10	6	CA
	66027327	4" Circular (17g)	5	4	CA
	66027328	3 1/2" x 1" Tubular (3g)	10	6	CA
	66027329	4 3/4" x 1 1/2" Tubular (7g)	5	4	CA


Gentle Adhesive


ALLEVYN Thin					
	Code	Description	Pcs/Pkg	Pkgs/Unit	Order Unit
	66047576	2" x 2 3/8"	10	10	CA
	66047578	4" x 4"	5	20	CA
	66047579	6" x 8"	3	16	CA


ALLEVYN Compression					
	Code	Description	Pcs/Pkg	Pkgs/Unit	Order Unit
	66047581	2" x 2 3/8"	10	10	CA
	66047583	4" x 4"	5	20	CA
	66047584	6" x 8"	3	16	CA


ALLEVYN Plus Cavity					
	Code	Description	Pcs/Pkg	Pkgs/Unit	Order Unit
	66047571	2" x 2 3/8"	10	10	CA
	66047573	4" x 4"	5	20	CA
	66047574	6" x 8"	3	16	CA


Secure Adhesive

ALLEVYN Island					
	Code	Description	Pcs/Pkg	Pkgs/Unit	Order Unit
	59764400	4" x 4"	5	10	CA
	59764500	5" x 5"	5	10	CA
	59764600	6" x 6"	5	10	CA
	59764700	6" x 5 3/8" Sacral	5	5	CA

ALLEVYN Adhesive					
	Code	Description	Pcs/Pkg	Pkgs/Unit	Order Unit
	66020043	3" x 3"	10	4	CA
	66020044	5" x 5"	10	4	CA
	66020045	7" x 7"	10	4	CA
	66000046	9" x 9"	10	8	CA

ALLEVYN Sacrum					
	Code	Description	Pcs/Pkg	Pkgs/Unit	Order Unit
	66000700	6 3/4" x 6 3/4"	10	8	CA
	66000451	9" x 9"	10	8	CA

ALLEVYN Plus Adhesive					
	Code	Description	Pcs/Pkg	Pkgs/Unit	Order Unit
	66000805	5" x 5"	10	5	CA
	66000806	7" x 7"	10	5	CA
	66000807	5" x 9"	10	5	CA

ALLEVYN Plus Adhesive Sacrum					
	Code	Description	Pcs/Pkg	Pkgs/Unit	Order Unit
	66000856	6 3/4" x 6 3/4"	10	5	CA
	66000857	9" x 9"	10	5	CA

Wound Management Customer Care Center www.allevyn.com
 Smith & Nephew, Inc. 1 800-876-1261 www.smith-nephew.com
 11775 Starkey Road T 727-392-1261
 P.O. Box 1970 F 727-392-6914
 Largo, Florida 33779
 USA